

North State Radiology

1720 Esplanade
Chico, CA 95926
(530) 898-0504
Fax: (530) 898-9647

AMICAS PACS Basic Workstation Viewer

Recommended System Specifications

August 2010

Hardware

- **CPUs** 2.66GHz Intel Dual Core CPU or similar
- **RAM** 2 GB
- **Video card** 128MB, True Color, 1280 x 1024 resolution or greater
- **CD Burner** Must have a functioning CD Burner (if Personal CD Burning feature will be used)
- **Internet Connection** DSL, Cable, or T1

Software

- **Operating System** Windows XP Professional, Windows Vista Business
- **Java Client** Java Version 6 Update 7, or higher
(will install automatically if user has local administrative privileges for the workstation)
- **Web Browser** Microsoft Internet Explorer version 7.0, or
Microsoft Internet Explorer v. 8.0 (compatibility view)
- **Anti-Virus Software** Full-featured (*i.e.*, “corporate” or “enterprise”) version of an industry-standard anti-virus program

Minimum System Specifications

- **CPUs** Pentium 4 or equivalent
- **RAM** 1GB
- **Operating System** Windows XP

***We also recommend the following security measures to ensure normal Java operation:**

--Use of Anti-Spyware and/or Personal Firewall software, with updated threat definitions.

--Installation of Critical and/or High Priority security updates via Windows Automatic Updates or the Microsoft Windows Update website.

Ata Alijani, M.D.

Frederick R. Brandon, D.O.

John C. Burnham, M.D.

Benjamin H. Cannon, M.D.

Donald E. Decker, M.D.

Gianluigi Fornari, M.D.

Garth L. Gregory, M.D.

Marcia J. Hinkle, M.D.

Robert P. Huberman, M.D.

Douglas M. Hughes, M.D.

Gregory J. Lauten, M.D.

Craig S. Letner, M.D.

J. Howard McIntosh, M.D.

Kevin M. Nash, M.D.

A. Elton Noyes, M.D.

Kenneth G. Petty, M.D.

George D. Rawley, M.D.

James F. Schlund, M.D.

Mauricio T. Schrader, M.D.

James F. Silverman, M.D.

Kiran K. Singh, M.D.

Gary A. Stromberg, M.D.

James F. Warren, M.D.